


# Disposable Electrosurgical Knife

## Indications

In endoscopic surgery, high-frequency current is used to excise and collect tissue samples in the digestive tract, and electrocoagulation to stop bleeding.



## Features

- ※ Gold coating on the blade effectively prevents sticking
- ※ I-shaped cutter head: quick dot marking, mucosal stripping
- ※ O-shaped cutter head: ceramic insulator design on the head, greatly improving surgical efficiency and cuttingsafety
- ※ L-shaped cutter head: rotatable cutter head, large area mucosa peeling
- ※ T-shaped cutter head: can be used for marking, hemostasis, pre-cutting, incision, peeling and other processes



## Specifications

Code	Cutter Model	Sheath O.D. ±0.2 ( mm )	Sheath Working Length ±100 (mm )	Cutter Extension Length ±0.5(mm)	Water Injection	Min. Access Channel Φ (mm)
All specifications	I	2.4	1600 1800 1950 2300	1.5/2/4	Yes/No	2.8
	O			2/4	Yes/No	
	T			2/4	Yes/No	
	L			4	No	